

# Reg3b Antibody, FITC conjugated



PACO34352

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Bactericidal C-type lectin which acts against several intestinal Gram-positive bacteria and Gram-negative bacteria. Lacks antibacterial activity against *S. typhimurium*. May play a role in protection against infection with *S. enteritidis* by inhibiting its translocation from the gut lumen into intestinal tissues and further extraintestinal tissues.

**Gene ID:**

Reg3b

**Uniprot**

P25031

**Synonyms:**

Regenerating islet-derived protein 3-beta (REG-3-beta) (Pancreatitis-associated protein 1) (Peptide 23) (REG-2) (Regenerating islet-derived protein III-beta) (Reg III-beta) [Cleaved into: Regenerating islet-derived protein 3-beta 16.5 kDa form; Regenerating islet-derived protein 3-beta 15 kDa form], Reg3b, Pap Pap1 Reg2

**Immunogen:**

Recombinant Rat Regenerating islet-derived protein 3-&beta; protein (27-175AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A