

# AAP1 Antibody



PACO34362

---

## Product Information

**Size:**

50ug

**Reactivity:**

Saccharomyces cerevisiae

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Positive effector of glycogen accumulation. May be involved in nutrient-sensing.

**Gene ID:**

AAP1

**Uniprot**

P37898

**Synonyms:**

Alanine/arginine aminopeptidase (EC 3.4.11. -), AAP1

**Immunogen:**

Recombinant Saccharomyces cerevisiae Alanine/arginine aminopeptidase protein (1-389AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A