

Major pollen allergen Car b 1 isoforms 1A and 1B Antibody

PACO34390



Product Information

Size:

50ug

Reactivity:

Carpinus betulus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:****Uniprot**

P38949

Synonyms:

Major pollen allergen Car b 1 isoforms 1A and 1B (Allergen Car b I) (allergen Car b 1)


Immunogen:

Recombinant Carpinus betulus Major pollen allergen Car b 1 isoforms 1A and 1B protein (2-160AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: Carpinus betulus (20µg). All lanes: Allergen Car b I antibody at 4µg/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 18 kDa. Observed band size: 18 kDa.