

Allergen Ara h 1, clone P17 Antibody



PACO34398

Product Information

Size:

50ug

Reactivity:

Peanut

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:****Uniprot**

P43237

Synonyms:

Allergen Ara h 1, clone P17 (Allergen Ara h I) (allergen Ara h 1)

Immunogen:

Recombinant Peanut Allergen Ara h 1, clone P17 protein (26-216AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A