

**Allergen Ara h 1, clone P17 Antibody, FITC
conjugated
PACO34400**



Product Information

Size:

50ug

Protein Background:

Reactivity:

Peanut

Gene ID:

Source:

Rabbit

Uniprot

P43237

Isotype:

IgG

Synonyms:

Allergen Ara h 1, clone P17 (Allergen Ara h I) (allergen Ara h 1)

Applications:

ELISA

Immunogen:

Recombinant Peanut Allergen Ara h 1, clone P17 protein (26-216AA).

Recommended dilutions:

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A