

Streptavidin Antibody



PACO34402

Product Information

Size:

50ug

Reactivity:

Streptomyces avidinii

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The biological function of streptavidin is not known. Forms a strong non-covalent specific complex with biotin (one molecule of biotin per subunit of streptavidin).

Gene ID:**Uniprot**

P22629

Synonyms:

Streptavidin

Immunogen:

Recombinant Streptomyces avidinii Streptavidin protein (25-183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A