

PACO34406

Product Information

Size:

50ug

Reactivity:

Cucurbita maxima

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May be involved in a redox system involving ascorbic acid, .

Gene ID:

AAO

Uniprot

P24792

Synonyms:

L-ascorbate oxidase (ASO) (Ascorbase) (EC 1.10.3.3), AAO


Immunogen:

Recombinant Cucurbita maxima L-ascorbate oxidase protein (31-579AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: Cucurbita maxima (40 μg/g, 20 μg/g, 10 μg/g, 5 μg/g). All lanes: AAO antibody at 3 μg/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 65 kDa. Observed band size: 65 kDa.