

# GV51\_0914 Antibody



PACO34414

---

## Product Information

**Size:**

50ug

**Reactivity:**

Gardnerella vaginalis 5-1

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

GV51\_0914

**Uniprot**

D6SZ78

**Synonyms:****Immunogen:**

Recombinant Gardnerella vaginalis 5-1 Modification methylase Accl protein (1-525AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A