

GV51_0914 Antibody, HRP conjugated



PACO34415

Product Information

Size:

50ug

Reactivity:

Gardnerella vaginalis 5-1

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GV51_0914

Uniprot

D6SZ78

Synonyms:

Immunogen:

Recombinant Gardnerella vaginalis 5-1 Modification methylase Accl protein (1-525AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A