

PACO34418

Product Information

Size:

50ug

Reactivity:

Hepatitis C virus genotype 1a

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P26664

Synonyms:

Genome polyprotein [Cleaved into: Core protein p21 (Capsid protein C) (p21); Core protein p19; Envelope glycoprotein E1 (gp32) (gp35); Envelope glycoprotein E2 (NS1) (gp68) (gp70); p7; Protease NS2-3 (p23) (EC 3.4.22. -); Serine protease NS3 (EC 3.4.21.98) (EC 3.6.1.15) (EC 3.6.4.13) (Hepacivirin) (NS3P) (p70); Non-structural protein 4A (NS4A) (p8); Non-structural protein 4B (NS4B) (p27); Non-structural protein 5A (NS5A) (p56); RNA-directed RNA polymerase (EC 2.7.7.48) (NS5B) (p68)]

Immunogen:

Recombinant Hepatitis C virus genotype 1a Genome polyprotein protein (192-325AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A