
Product Information

Size:

50ug

Reactivity:

Actinoplanes utahensis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the hydrolysis of the palmitoyl moiety of the antifungal antibiotic, aculeacin-A, giving a hexapeptide moiety and a long chain fatty acid, .

Gene ID:

aac

Uniprot

P29958

Synonyms:

Aculeacin-A acylase (EC 3.5.1. -) [Cleaved into: Aculeacin-A acylase small subunit; Aculeacin-A acylase large subunit], aac

Immunogen:

Recombinant Actinoplanes utahensis Aculeacin-A acylase protein (35-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A