

R107.2 Antibody, HRP conjugated



PACO34431

Product Information

Size:

50ug

Reactivity:

Caenorhabditis elegans

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the dehydration of the S-form of NAD(P)HX at the expense of ATP, which is converted to ADP. Together with NAD(P)HX epimerase, which catalyzes the epimerization of the S- and R-forms, the enzyme allows the repair of both epimers of NAD(P)HX, a damaged form of NAD(P)H that is a result of enzymatic or heat-dependent hydration.

Gene ID:

R107.2

Uniprot

P32740

Synonyms:

ATP-dependent (S)-NAD(P)H-hydrate dehydratase (EC 4.2.1.93) (ATP-dependent NAD(P)HX dehydratase)

Immunogen:

Recombinant Caenorhabditis elegans ATP-dependent (S)-NAD(P)H-hydrate dehydratase protein (1-307AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A