

# E6 Antibody



PACO34438

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human papillomavirus type 52

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcriptional transactivator. Binds double-stranded DNA.

**Gene ID:**

E6

**Uniprot**

P36814

**Synonyms:**

Protein E6, E6

**Immunogen:**

Recombinant Human papillomavirus type 52 Protein E6 protein (1-148AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A