

CYCB1-1 Antibody




PACO34482

Product Information

Size: 50ug	Protein Background:
Reactivity: Arabidopsis thaliana	Gene ID: CYCB1-1
Source: Rabbit	Uniprot P30183
Isotype: IgG	Synonyms: Cyclin-B1-1 (Cyc1-At) (G2/mitotic-specific cyclin-B1-1) (CycB1;1), CYCB1-1, CYC1
Applications: ELISA, WB	Immunogen: Recombinant Arabidopsis thaliana Cyclin-B1-1 protein (1-428AA).
Recommended dilutions: ELISA:1:2000-1:10000, WB:1:1000-1:5000	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: Arabidopsis thaliana (40µg, 20µg). All lanes: CYCB1-1 antibody at 1.8µg/ml. Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 49 kDa. Observed band size: 49 kDa.