

# Ribosome-inactivating protein bryodin I Antibody PACO34486



---

## Product Information

**Size:**

50ug

**Reactivity:**

Bryonia dioica

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Ribosome-inactivating protein of type 1, inhibits protein synthesis in animal cells.

**Gene ID:****Uniprot**

P33185

**Synonyms:**

Ribosome-inactivating protein bryodin I (EC 3.2.2.22) (BD1) (rRNA N-glycosidase)

**Immunogen:**

Recombinant Bryonia dioica Ribosome-inactivating protein bryodin I protein (24-247AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A