

PACO34498

Product Information

Size:

50ug

Reactivity:

Bacillus licheniformis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Hydrolyzes the cell wall of *M. luteus* more efficiently than that of *B. licheniformis* and *B. subtilis*. The C-terminal region, including the repeats, determines substrate specificity.

Gene ID:

cwIM

Uniprot

P37134

Synonyms:

N-acetylmuramoyl-L-alanine amidase CwIM (EC 3.5.1.28) (Autolysin) (Cell wall hydrolase), cwIM

Immunogen:

Recombinant *Bacillus licheniformis* N-acetylmuramoyl-L-alanine amidase CwIM protein (1-253AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A