

# PMAP36 Antibody



PACO34542

---

## Product Information

**Size:**

50ug

**Reactivity:**

Pig

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Exerts antimicrobial activity against both Gram-positive and negative bacteria. Its activity appears to be mediated by its ability to damage bacterial membranes.

**Gene ID:**

PMAP36

**Uniprot**

P49931

**Synonyms:**

Antibacterial peptide PMAP-36 (Myeloid antibacterial peptide 36), PMAP36

**Immunogen:**

Recombinant Pig Antibacterial peptide PMAP-36 protein (130-166AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A