

PACO34550

Product Information

Size:

50ug

Reactivity:

Helicobacter pylori

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be necessary for the transcription, folding, export, or function of the cytotoxin.

Gene ID:

cagA

Uniprot

P55980

Synonyms:

Cytotoxicity-associated immunodominant antigen (120 kDa protein) (CAG pathogenicity island protein 26), cagA, cag26 cai

Immunogen:

Recombinant Helicobacter pylori Cytotoxicity-associated immunodominant antigen protein (1-294AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A