

# Selenot Antibody



PACO34586

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Selenot

**Uniprot**

P62342

**Synonyms:**

Thioredoxin reductase-like selenoprotein T (SelT) (EC 1.8.1.9), Selenot, Selt

**Immunogen:**

Recombinant Mouse Thioredoxin reductase-like selenoprotein T protein (20-195AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A