

PACO34590

Product Information

Size:

50ug

Reactivity:

Zaire ebolavirus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

sGP seems to possess an anti-inflammatory activity as it can reverse the barrier-decreasing effects of TNF alpha. Might therefore contribute to the lack of inflammatory reaction seen during infection in spite the of extensive necrosis and massive virus production. Does not seem to be involved in activation of primary macrophages. Does not seem to interact specifically with neutrophils.

Gene ID:

GP

Uniprot

P60171

Synonyms:

Pre-small/secreted glycoprotein (pre-sGP) [Cleaved into: Small/secreted glycoprotein (sGP); Delta-peptide], GP

Immunogen:

Recombinant Zaire ebolavirus Pre-small/secreted glycoprotein protein (33-364AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A