

H2-K1 Antibody, HRP conjugated



PACO34627

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

H2-K1

Uniprot

P01902

Synonyms:

H-2 class I histocompatibility antigen, K-D alpha chain (H-2K(D)), H2-K1, H2-K

Immunogen:

Recombinant Mouse H-2 class I histocompatibility antigen, K-D α chain protein (22-305AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A