

Minor capsid protein VP2 Antibody



PACO34634

Product Information

Size:

50ug

Reactivity:

JC polyomavirus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P03095

Synonyms:

Minor capsid protein VP2 (Minor structural protein VP2)

Immunogen:

Recombinant JC polyomavirus Minor capsid protein VP2 protein (2-344AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A