

Histidine-rich protein PFHRP-II Antibody



PACO34686

Product Information

Size:

50ug

Reactivity:

Plasmodium falciparum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P05227

Synonyms:

Histidine-rich protein PFHRP-II

Immunogen:

Recombinant Plasmodium falciparum Histidine-rich protein PFHRP-II protein (31-289AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A