

L2 Antibody



PACO34694

Product Information

Size:

50ug

Reactivity:

Human papillomavirus type 18

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

L2

Uniprot

P06793

Synonyms:

Minor capsid protein L2, L2

Immunogen:

Recombinant Human papillomavirus type 18 Minor capsid protein L2 protein (1-462AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A