

Fcgr3 Antibody



PACO34718

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for the Fc region of complexed immunoglobulins gamma. Low affinity receptor.

Gene ID:

Fcgr3

Uniprot

P08508

Synonyms:

Low affinity immunoglobulin gamma Fc region receptor III (IgG Fc receptor III) (Fc-gamma RIII) (FcRIII) (CD antigen CD16), Fcgr3

Immunogen:

Recombinant Mouse Low affinity immunoglobulin γ Fc region receptor III protein (31-215AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A