

PACO34722

---

## Product Information

**Size:**

50ug

**Reactivity:**

Enterobacteria phage 933W

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

The B subunit is responsible for the binding of the holotoxin to specific receptors on the target cell surface, such as globotriaosylceramide (Gb3) in human intestinal microvilli.

**Gene ID:**

stxB2

**Uniprot**

P09386

**Synonyms:**

Shiga-like toxin 2 subunit B (SLT-2 B subunit) (SLT-2b) (SLT-IIb) (Verocytotoxin 2 subunit B) (Verotoxin 2 subunit B), stxB2, stx2B

**Immunogen:**

Recombinant Enterobacteria phage 933W Shiga-like toxin 2 subunit B protein (20-89AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A