

Ribosome-inactivating protein alpha-trichosanthin Antibody

PACO34726



Product Information

Size:

50ug

Reactivity:

Trichosanthes kirilowii

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inactivates eukaryotic 60S ribosomal subunits.

Gene ID:**Uniprot**

P09989

Synonyms:

Ribosome-inactivating protein alpha-trichosanthin (Alpha-TCS) (EC 3.2.2.22) (rRNA N-glycosidase)

Immunogen:

Recombinant Trichosanthes kirilowii Ribosome-inactivating protein alpha-trichosanthin protein (24-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A