

mpt64 Antibody, HRP conjugated



PACO34759

Product Information

Size:

50ug

Reactivity:

Mycobacterium tuberculosis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

mpt64

Uniprot

P9WIN9

Synonyms:

Immunogenic protein MPT64 (Antigen MPT64), mpt64

Immunogen:

Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A