

# mpt64 Antibody, FITC conjugated



PACO34760

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mycobacterium tuberculosis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

mpt64

**Uniprot**

P9WIN9

**Synonyms:**

Immunogenic protein MPT64 (Antigen MPT64), mpt64

**Immunogen:**

Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A