

PACO34770

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

The RecF protein is involved in DNA metabolism; it is required for DNA replication and normal SOS inducibility. RecF binds preferentially to single-stranded, linear DNA. It also seems to bind ATP.

Gene ID:

recF

Uniprot

P0A7H0

Synonyms:

DNA replication and repair protein RecF, recF, uvrF

Immunogen:

Recombinant Escherichia coli DNA replication and repair protein RecF protein (2-357AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A