

PACO34794

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

crp

Uniprot

P0ACJ8

Synonyms:

cAMP-activated global transcriptional regulator CRP (Catabolite activator protein) (CAP) (Catabolite gene activator) (cAMP receptor protein) (CRP) (cAMP regulatory protein), crp, cap csm

Immunogen:

Recombinant Escherichia coli cAMP-activated global transcriptional regulator CRP protein (1-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A