

PACO34806

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

integral component of membrane, integral component of plasma membrane, plasma membrane, endopeptidase activity, peptidase activity, serine-type endopeptidase activity, toxic substance binding, protein processing, proteolysis, signal peptide processing.

Gene ID:

lepB

Uniprot

P00803

Synonyms:

Signal peptidase I (SPase I) (EC 3.4.21.89) (Leader peptidase I), lepB

Immunogen:

Recombinant Escherichia coli Signal peptidase I protein (78-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A