

Hbb-b2 Antibody



PACO34810

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Hbb-b2

Uniprot

P02089

Synonyms:

Hemoglobin subunit beta-2 (Beta-2-globin) (Hemoglobin beta-2 chain) (Hemoglobin beta-minor chain), Hbb-b2

Immunogen:

Recombinant Mouse Hemoglobin subunit beta-2 protein (2-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A