

# Hbb-b2 Antibody, Biotin conjugated



PACO34813

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Hbb-b2

**Uniprot**

P02089

**Synonyms:**

Hemoglobin subunit beta-2 (Beta-2-globin) (Hemoglobin beta-2 chain) (Hemoglobin beta-minor chain), Hbb-b2

**Immunogen:**

Recombinant Mouse Hemoglobin subunit beta-2 protein (2-147AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A