

# BLLF1 Antibody



PACO34826

---

## Product Information

**Size:**

50ug

**Reactivity:**

Epstein-Barr virus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Initiates virion attachment to host B-lymphocyte cell, leading to virus entry. Acts by binding to host CR2 at the surface of B-lymphocytes, facilitating the binding of viral glycoprotein gp42 to HLA class II molecules. Attachment triggers virion-host membrane fusion and invasion of the host cell.

**Gene ID:**

BLLF1

**Uniprot**

P03200

**Synonyms:**

Envelope glycoprotein GP350 (Membrane antigen) (MA)

**Immunogen:**

Recombinant Epstein-Barr virus Envelope glycoprotein GP350 protein (1-493AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A