

# Capsid protein VP1 Antibody



PACO34862

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human parvovirus B19

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

**Uniprot**

P07299

**Synonyms:**

Capsid protein VP1 (Coat protein VP1) [Cleaved into: Capsid protein VP2 (Coat protein VP2)]

**Immunogen:**

Recombinant Human parvovirus B19 Capsid protein VP2 protein (228-781AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A