

Capsid protein VP1 Antibody, Biotin conjugated



PACO34865

Product Information

Size:

50ug

Reactivity:

Human parvovirus B19

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:****Uniprot**

P07299

Synonyms:

Capsid protein VP1 (Coat protein VP1) [Cleaved into: Capsid protein VP2 (Coat protein VP2)]

Immunogen:

Recombinant Human parvovirus B19 Capsid protein VP2 protein (228-781AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A