

PACO34866

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Apo(a) is the main constituent of lipoprotein(a) (Lp(a)). It has serine proteinase activity and is able of autoproteolysis. Inhibits tissue-type plasminogen activator 1. Lp(a) may be a ligand for megalin/Gp 330.

Gene ID:

LPA

Uniprot

P08519

Synonyms:

Apolipoprotein(a) (Apo(a)) (Lp(a)) (EC 3.4.21. -), LPA

Immunogen:

Recombinant Human Apolipoprotein(a) protein (3999-4112AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A