

Trypsin inhibitor DE-3 Antibody, HRP conjugated



PACO34879

Product Information

Size:

50ug

Reactivity:

Erythrina caffra

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibition of trypsin.

Gene ID:

Uniprot

P09943

Synonyms:

Trypsin inhibitor DE-3

Immunogen:

Recombinant Erythrina caffra Trypsin inhibitor DE-3 protein (1-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A