

inhA Antibody



PACO34902

Product Information

Size:

50ug

Reactivity:

Mycobacterium tuberculosis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

inhA

Uniprot

P9WGR1

Synonyms:

Enoyl-[acyl-carrier-protein] reductase [NADH] (ENR) (Enoyl-ACP reductase) (EC 1.3.1.9) (FAS-II enoyl-ACP reductase) (NADH-dependent 2-trans-enoyl-ACP reductase), inhA

Immunogen:

Recombinant Mycobacterium tuberculosis Enoyl-[acyl-carrier-protein] reductase [NADH] protein (1-269AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A