

PACO34918

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in DNA repair and RecF pathway recombination.

Gene ID:

recO

Uniprot

P0A7H3

Synonyms:

DNA repair protein RecO (Recombination protein O), recO

Immunogen:

Recombinant Escherichia coli DNA repair protein RecO protein (1-242AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A