

entC1 Antibody



PACO34938

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Staphylococcal enterotoxins cause the intoxication staphylococcal food poisoning syndrome. The illness characterized by high fever, hypotension, diarrhea, shock, and in some cases death.

Gene ID:

entC1

Uniprot

P01553

Synonyms:

Enterotoxin type C-1 (SEC1), entC1

Immunogen:

Recombinant Staphylococcus aureus Enterotoxin type C-1 protein (28-266AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A