

H2-D1 Antibody



PACO34946

Product Information

Size:

50ug

Reactivity:

Mouse, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Involved in the presentation of foreign antigens to the immune system.

Gene ID:

H2-D1

Uniprot

P01900

Synonyms:

H-2 class I histocompatibility antigen, D-D alpha chain (H-2D(D)), H2-D1


Immunogen:

Recombinant Mouse H-2 class I histocompatibility antigen, D-D α chain protein (25-311AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: H2-D1 antibody at 2 μ g/ml + A549 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 42 kDa. Observed band size: 42 kDa.