

# Parvalbumin beta Antibody, HRP conjugated



PACO34951

---

## Product Information

**Size:**

50ug

**Reactivity:**

Gadus morhua subsp

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

In muscle, parvalbumin is thought to be involved in relaxation after contraction. It binds two calcium ions.

**Gene ID:**

**Uniprot**

P02622

**Synonyms:**

Parvalbumin beta (Allergen Gad c I) (Allergen M) (allergen Gad c 1)

**Immunogen:**

Recombinant Gadus morhua subsp Parvalbumin &beta; protein (1-113AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A