

E6 Antibody



PACO34958

Product Information

Size:

50ug

Reactivity:

Human papillomavirus type 16

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

E6

Uniprot

P03126

Synonyms:

Protein E6, E6

Immunogen:

Recombinant Human papillomavirus type 16 Protein E6 protein (1-158AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A