

# d2hgdh Antibody



PACO34982

---

## Product Information

**Size:**

50ug

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the oxidation of D-2-hydroxyglutarate to alpha-ketoglutarate.

**Gene ID:**

d2hgdh

**Uniprot**

A1L258

**Synonyms:**

D-2-hydroxyglutarate dehydrogenase, mitochondrial (EC 1.1.99. -), d2hgdh

**Immunogen:**

Recombinant Zebrafish D-2-hydroxyglutarate dehydrogenase, mitochondrial protein (56-533AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A