

amiF Antibody, HRP conjugated



PACO35011

Product Information

Size:

50ug

Reactivity:

Bradyrhizobium sp.

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Is an aliphatic amidase with a restricted substrate specificity, as it only hydrolyzes formamide.

Gene ID:

amiF

Uniprot

A4Z3G9

Synonyms:

Formamidase (EC 3.5.1.49) (Formamide amidohydrolase), amiF

Immunogen:

Recombinant Bradyrhizobium sp. Formamidase protein (1-337AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A