

PACO35054

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

blood microparticle, external side of plasma membrane, extracellular exosome, immunoglobulin complex, circulating, antigen binding, immunoglobulin receptor binding, B cell receptor signaling pathway, complement activation, classical pathway, defense response to bacterium, Fc-epsilon receptor signaling pathway.

**Gene ID:**

IGLL5

**Uniprot**

B9A064

**Synonyms:**

Immunoglobulin lambda-like polypeptide 5 (G lambda-1) (Germline immunoglobulin lambda 1), IGLL5

**Immunogen:**

Recombinant Human Immunoglobulin lambda-like polypeptide 5 protein (36-214AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A