

**Phospholipase D LamSicTox-alphaIC1 Antibody,
FITC conjugated**
PACO35060



Product Information

Size:

50ug

Reactivity:

Loxosceles amazonica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the hydrolysis of sphingomyelin. May also acts on other phosphatidyl esters. Induces complement-dependent hemolysis, dermonecrosis, blood vessel permeability and platelet aggregation.

Gene ID:

Uniprot

C0JAZ9

Synonyms:

Phospholipase D LamSicTox-alphaIC1 (PLD) (EC 3.1.4.4) (Dermonecrotic toxin) (Sphingomyelin phosphodiesterase D) (SMD) (SMase D) (Sphingomyelinase D) (Fragment)

Immunogen:

Recombinant Loxosceles amazonica Phospholipase D LamSicTox-alphaIC1 protein (1-273AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A