

ipaH9.8 Antibody



PACO35078

Product Information

Size:

50ug

Reactivity:

Shigella flexneri serotype X

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Effector proteins function to alter host cell physiology and promote bacterial survival in host tissues. This protein is an E3 ubiquitin ligase that interferes with host's ubiquitination pathway and modulates the acute inflammatory responses, thus facilitating bacterial colonization within the host cell.

Gene ID:

ipaH9.8

Uniprot

D2AJU0

Synonyms:

E3 ubiquitin-protein ligase ipaH9.8 (EC 2.3.2.27) (Invasion plasmid antigen ipaH9.8) (RING-type E3 ubiquitin transferase ipaH9.8), ipaH9.8

Immunogen:

Recombinant Shigella flexneri serotype X E3 ubiquitin-protein ligase ipaH9.8 protein (1-545AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A